Easy Screw Removal KIT

User Manual

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English Easy Screw Removal KIT User Manual

1. Description

ESR(Easy Screw Removal) KIT is a kit with functions for abutment & abutment screw fracture removal, slipped connection part removal, and implant's internal screw thread restoration. The ESR KIT consists of abutment & abutment screw fracture removal device, slipped connection part removal device, and implant's internal screw thread-restoring device.

2. Purpose of Product

The purpose of this product is to remove the abutment & abutment screw fracture, remove the slipped connection part, and restore the implant's internal screw thread.

3. How to Use**• A guide to removing damaged screws**

1. Prepare a guide and reverse drill suitable for the fixture.
 - For guide and reverse drill, choose from either the Short (option) or Long type, considering the jawbones of the patients.
2. Attach the guide to the fixture and fix the guide with ESR handle.
 - If the operator attaches the guide window so that it is fully visible to the operator, the operation process can be seen and water can be injected inside.
 - For the external type, use heavy impression material to fix the guide firmly.
3. Set the engine for operation at 50rpm and reverse rotation mode, and attach the reverse drill to the hand-piece.
4. Rotate the reverse drill slowly at 30~50rpm in reverse direction and the damaged screw will back out for removal.
 - ※ If the damaged screw is not removed with the reverse drill, as in procedure no. 4,
5. Prepare the screw removal drill (hereinafter referred to as 'SR drill') suitable to the fixture.
 - For the SR drill, select the Short (option) or Long type, considering the jawbones of the patients.
6. Set the engine for operation at 1,500rpm and reverse rotation mode and attach the SR drill to the hand-piece.
7. Remove the damaged screw while rotating the SR drill in reverse direction and injecting water.
8. When drilling with the force of 5~10N (under 1kg), pump gently and drill until the red band of the sack disappears.

9. Insert the screw removal tip inside of guide and rotate it slowly in the reverse direction to remove the damaged screw.
10. If the internal screw thread of fixture is damaged, use Re-tap to restore the screw thread.
11. Inject water sufficiently inside the fixture to remove foreign materials such as remaining metals by suction.

• Fractured Abutment Removal Guide

1. In case of 1-Piece Abutment fracture, a hole is formed using dental bur and connected by reverse rotation using the Abutment removal Tip; the abutment is then removed by shaking using forceps.
2. In case of 2-Piece Abutment fracture, it is connected to the internal hole of the abutment by reverse rotation using the Abutment removal Tip; the abutment is then removed by shaking using forceps.

• Abutment Screw Hex Slip Removal Guide

1. In case of slip in the 1.2 Screw Hex part, the Hex part is slightly grinded using a small bur, and then connected using the Abutment Removal Tip for the Mini option and reverse-rotated for the removal.

• Cover Screw Hex Slip Removal Guide

1. In case of slip in the Cover Screw Hex connection part, it can be removed with a slot driver after making a straight “—” groove using a dental Ø 0.8 bur.

4. Notices when using

- Please read how to use completely before the operation.
 - This product is a medical device. Make sure to process moist heat sterilization in autoclave before using (temperature : 132°C, duration : 15 minutes).
1. Screw removal drill (SR drill) is limited for one-time use and discarded after the use. Re-using is prohibited.
 2. SR drill should be attached to the guide to prevent the damage of screw thread inside the fixture and drill blade, due to shaking.
 3. The recommended speed of SR Drill is 1,200~1,500rpm. Use low-or high-speed hand-piece and do not use at more than 2,000rpm.
 4. As metal powder may be created in the process of removing the damaged screw, inject water sufficiently after the operation and remove the foreign

English Easy Screw Removal KIT User Manual

- materials by suction.
5. The recommended speed of the reverse drill is 30~50rpm and the use is limited to 10-times.

5. Maintenance and Care

- Separate all used tools immediately after surgery, after cleaning thoroughly in the running water using a cleaning brush and drying. Store at room temperature.
- Do not leave the instruments in a place vulnerable to contamination.
- Do not use hydrogen peroxide as a disinfectant or cleaner (it may cause damage or discoloration of the TiN Coating, laser marking and color coding).
- This product comes with a 1-year warranty from the date of purchase.
- Because Screw removal Tip and Screw removal Drill are one-time use only, they should be discarded after use and re-using is prohibited.

6. Device information

For device information such as part name, part number and specifications, manufacturing date, lot number, and expiration date refer to the label on the package.

* **Note** : ※ For additional details, refer to catalogs and operation manuals. For general information on implants, refer to the catalog, prosthetic manual, surgical manual and other related data.

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Rx only For USA Only : Federal law restricts this device to sale by or on the order of a dentistDEUTSCHE OSTEM GmbH