OSSTEM[®]
IMPLANT

Lateral Approach - Sinus KIT (LAS-KIT)

User Manual

OSSTEM[®]
IMPLANT

Storage Condition :
Dry place at room temperature(1°C~30°C)

Caution, Consult
accompanying
Documents

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<http://www.osstem.com>

Rx only

For USA Only : Federal law restricts this
device to sale by or on the order of a dentist

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English

■ Product Introduction

Lateral Approach–Sinus KIT(LAS-KIT) is an operation tool optimized for the Lateral Approach which can elevate the floor of maxillary sinus membrane fast and safely. Please refer to our product catalog or brochure for detail information.

■ Washing KIT

After performing a dental surgery, wash the tools using the following procedures;

1. During surgery, keep the used tools in a saline solution or distilled water.
2. When the surgery is completed, soak all the tools used in the surgery in alcohol for washing. Caution: Do not wash the tools with hydrogen peroxide. Exposure to hydrogen peroxide may result in the discoloration of the laser marking or anodizing.
3. Wash the tools with distilled or clean flowing water to remove blood stain and other foreign matter.
4. Remove the moisture from the tools and set them in the KIT case according to the color coding for easy and correct set-up.
5. After setting-up, dry up the KIT in an autoclave at 132°C for at least 15 minutes, and then keep it at room temperature.

■ Cautions for Use

1. Select the appropriate drill type and drill diameter for the window formation according to the surgical plan.
2. When using the Dome drill or the Core drill, the initial drill depth should be less than the lateral wall thickness on the panorama or CT by 0.5 mm. Use the stopper that is connected to the drill. The side wall drill can be connected with the CAS-KIT Stopper if you need.

3. When the drilling, do not apply excessive force to the stopper when is has reached the desired position. Do not over drilling toward the sinus membrane.
4. When using the Dome drill or the Side Wall drill, avoid excessive lateral force.
5. The recommended cutting speed of the Dome drill, Core drill, and Side Wall drill is 1,200 - 1,500 rpm.
6. After making a window formation, membrane elevation can be performed by using the Sinus Lift Kit.
7. Keep the tools away from a location where there is a possibility of infection or contamination.
8. The products are medical devices; therefore, they must be sterilized in an autoclave at 132°C for at least 15 minutes before use.
9. Do not use hydrogen peroxide for sterilization or cleaning, otherwise, the laser marking, anodizing, or color coding can be damaged or discolored.
10. The recommended usage times of the Dome drill, Core drill, Side Wall drill, and Stopper is 50 times. However, this may vary depending on the conditions of use.
11. The warranty period of this product is one (1) year from the date of purchase.

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