

Periosteal Elevator (24G)

- Lifting mucosal periosteum after gingival tissue incision
- W : 4.2/4.0mm

EP24G-W

Minesota Retractor

- Securing a clear view by pilling the mouth, cheeks, etc.

RTCRM-W

Extension Hose Adapter

- Adapter for chair suction connection

SNKHA-W

Periosteal Elevator (Selden)

- Lifting mucosal periosteum after gingival tissue incision
- W : 10/13mm

EP23-W

Needle Holder (Crile-Wood, TC)

- Straight
- Tungsten carbide-treated beak
- L : 150mm (±5)

NHC150TC-W

Periodontal Chisel

- Bone removal and formation
- W : 5.0mm

CHCO2-W

Tissue Forcep (ADSON)

- Used to hold soft tissue
- No projections inside the beak
- L : 120mm (±5)

PT41-W

Extension Hose

- Extension hose for chair suction connection
- Autoclave can be used
- Transparent silicone material

SNKHS-1-W

Dr.Cho's Instrument Pouch

- Used for storing and sterilization of instruments
- L : 550 X 400mm

WPB-W

Titanium Suction Tip

- D(Inner Diameter) : 3.0mm

SN3TI-W

